

*C. Pena¹, D. Anderson¹, A. Apresyan¹, A. Bornheim¹, J. Duarte¹, A. Ronzhin², M. Spiropulu¹, S. Xie¹.
California Institute of Technology¹
Fermilab²*

Precision Timing Calorimeter for High Energy Physics

SCINT 2015



Cristián Peña
Caltech



Outline

- ❖ High Energy Physics motivation for pursuing precision timing in calorimeters
- ❖ Physical processes that affect precision timing (TOF resolution)
- ❖ Experimental beam test results (FTBF, CERN)
 - ❖ LYSO-based sampling calorimeter
 - ❖ LYSO-Tungsten Shashlik calorimeter
- ❖ Laser tests at Caltech
- ❖ Summary



Cristián Peña
Caltech



Physics Motivation for Precision Timing

High-Luminosity Large Hadron Collider (HL-LHC)

collect more data (3 ab^{-1}) \Rightarrow higher instantaneous luminosity \Rightarrow higher average primary interaction per bunch crossing (PU)



Goals

Precise measurements of Higgs properties:


- couplings, tensor structure, rare decay
- role in EW-SB ($W_L W_L$ scattering)

Explore signals of new physics:

- Supersymmetry, DM direct production
- Many analyses require jets and MET



Cristián Peña
Caltech





Physics Motivation for Precision Timing

Challenges

Collect more data to increase LHC reach:

- Pileup interactions up to 140
- Key measurements will be affected by this harsh environment



Cristián Peña
Caltech




Precision Timing Calorimeters

Calorimeter with time resolution of the $O(20-30)$ ps

- Allows a $H \rightarrow \gamma\gamma$ vertex reconstruction with approximately 1 cm resolution
- Reduces pileup (PU) energy by a factor of 5-10

Possible physics applications of timing information

- *Object level*: identify forward PU jets (improve VBF higgs and WW scattering)
- *Single hit level*: e.g timing-base ECAL clustering cleaning
- Spatially separate overlapping vertices that corresponds to different time



Cristián Peña
Caltech



Precision Timing Using Crystal


TOF resolution driven by a number of different effects

- Main effects can be approximately factorize
- EM shower development (t_I): shown to be less than 13 ps
(A. Ronzhin et. al. NIM-A, vol 749 p.65-73, new paper accepted by NIM-A)

We focus our studies on scintillation and transit time (t_S and t_T)

Setup allows to control:

Photodetector jitter (t_P) at the 10 ps level and DAQ resolution (t_D) to about 6ps




Cristián Peña
Caltech



Scintillation Light Time Spectrum

- Scintillating crystals are often classified according to its time decay constant (τ_d). This is of the order of tens of ns; PWO, LYSO ≈ 40 ns
- In a calorimeter: timing information is obtained from the leading edge of signal; thus the **rise time of the light output** is the key quantity
- LYSO scintillation light properties:
 - Light output rise time $t_R < 75$ ps - Light yield $\approx 35,000$ photons/MeV


See: S Seifert, J.H.L Steenbergen, H.T van Dam and D.R Schaart, 2012 *JINST* 7 P09004. doi:10.1088/1748-0221/7/09/P09004



Cristián Peña
Caltech




LYSO-based Sampling Calorimeter



All measurements were carried out using electron beams

All measurement carried out at the Fermilab Test Beam Facility (FTBF)

DAQ = DRS4 wave digitizer + laptop



Goal: study the effect of scintillation light by reducing optical transport.
Electromagnetic Shower developed in the lead radiator, then sampled by the LYSO crystal





$$\text{TOF} = t_1 - t_0 \quad (1)$$



Cristián Peña
Caltech



LYSO-based Sampling Calorimeter Results



- Measure TOF resolution for different electron energies (4, 8, 16, 32 GeV)
- Measure 33 ps TOF resolution for 32 GeV electrons
- TOF resolution: width of gaussian fit to the TOF distribution (σ)




Cristián Peña
Caltech

NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



LYSO-based Sampling Calorimeter Results II



- Observe a $1/\sqrt{E}$ behavior for TOF resolution
- Fit the TOF resolution distribution to the superposition of a stochastic term and a constant term
- **Summary: 20-30 ps TOF resolution goal is achievable for ~ 50 GeV e/γ objects; provided enough light collection**




Cristián Peña
Caltech



NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



Increased Light Path Complexity





- R&D for HL-LHC on LYSO/W shashlik calorimeter
- Radiation hard in HL-LHC up to 3 ab^{-1}
- Energy resolution of $10\%/\sqrt{E} \oplus 1\%$
- Energy resolution performance shown in test beam on a 4x4 matrix
- Use single shashlik cell to test timing performance



Cristián Peña
Caltech



LYSO-Tungsten Shashlik Calorimeter





Goal: study the effect of optical transit by increasing the distance traveled by the light
Goal: measure the TOF resolution using WLS fiber and side readout (a la sampling calorimeter)
Build single shashlik module with WLS fibers and MCP-PMT readout



Cristián Peña
Caltech



Shashlik WLS fibers Readout Results I



Test the response of different WLS fiber, particularly the rise time. Measure rise times to be ~ 7.1 ns and ~ 2.4 ns for the Y11 (kuraray) and DSB1 WLS fibers respectively.

In terms of TOF resolution the faster (DSB1) fiber performs better.




Cristián Peña
Caltech

NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



Shashlik Timing Performance



- Performance of LYSO based calorimeters scales with rise time
- Main difference is due to WLS
- **Summary:** few tens of ps TOF resolution is achievable using LYSO based calorimeters for e/γ objects of ~ 100 GeV




Cristián Peña
Caltech

NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



Current Systematic Limits

- Extended measurement up to 150 GeV suggests that systematic limit is of the order of 10 ps
- Fit yields constant term compatible with reference time resolution of 15 ps



Cristián Peña
Caltech





Light Sensing Performance

- Differential time resolution of two MCP-PMTs
- DRS4 read out
- Illuminated with a 50 ps FWHM laser
- $\sigma_t = 7.2$ ps
- Multi-photon timing performance of MCP-PMT approaches read-out limit

Hamamatsu MCP-PMT

Rise Time [Ⓞ]	150	—	ps
Fall Time [Ⓜ]	360	—	ps
I.R.F. (FWHM) [Ⓛ]	45 [Ⓜ]	—	ps
T.T.S. (FWHM)	—	25 [Ⓞ]	ps



Cristián Peña
Caltech



Summary

- ❖ Carried out measurement in beam test with standalone LYSO-based sampling and LYSO-Tungsten Shashlik calorimeter prototypes.
- ❖ The 20-30 ps TOF resolution goal is achieved using a LYSO-based sampling calorimeter with MCP-PMT and DRS4 read out.
- ❖ Study the optical transit time effect on TOF resolution by using different WLS fiber readout and by direct optical coupling to sides of the Shashlik cell.
- ❖ Measure $\sim 50/100$ ps TOF resolution for a single Shashlik cell using side/WLS(DBS1) read out
- ❖ A large scale calorimeter with a time resolution of few tens of ps for pileup mitigation seems achievable




Cristián Peña
Caltech





Backups



Cristián Peña
Caltech



Shashlik WLS fibers Readout Results II



Use same TOF algorithm as in the LYSO-sampling calorimeter.

Measure a $\sim 4.5\%$ energy resolution using a single shashlik cell readout by WLS fiber and MCP-PMTs

Observe a $1/\sqrt{E}$ behavior for the TOF resolution, measure ~ 100 ps time resolution for 32 GeV e^-





Cristián Peña
Caltech

NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



Shashlik Side Readout



Alternative timing readout:
direct coupling of
photodevice to the edges of
the shashlik cell. Decreases
transit time; Exposing two
LYSO tiles

Similar approach as a sampling calorimeter by direct optical coupling of MCP-PMTs at $\sim 5X_0$.
Reduces transit time jitter with at the expense of collecting less light.
Energy measurement obtained with WLS fibers





Cristián Peña
Caltech

NIM-A, vol 794 p. 7-4.
doi:10.1016/j.nima.2015.04.013



Shashlik Side Readout Results



Measure a faster rise time from direct side coupling readout compared to that of the WLS fibers.

Observe a $1/\sqrt{E}$ behavior for the TOF resolution, measure ~ 50 ps time resolution for 32 GeV e^-



Cristián Peña
Caltech



Hamamatsu MCP-PMT R3809U-50

http://www.hamamatsu.com/resources/pdf/etd/R3809U-50_TPMH1067E10.pdf



Hamamatsu MCP-PMT

Rise Time [Ⓞ]	150	—	ps
Fall Time [Ⓢ]	360	—	ps
I.R.F. (FWHM) ^①	45 ^②	—	ps
T.T.S. (FWHM)	—	25 [Ⓞ]	ps

Photek 240

http://www.photek.com/pdf/datasheets/detectors/DS006_Photomultipliers.pdf



Cristián Peña

Caltech



Ability to survive up to 3000 fb^{-1} , and to 2035

- predicted neutron fluence of about 10^{16} n/cm^2 in forward regions (CMS)
- ionizing dose up to 150 Mrad in CMS electromagnetic calorimeter ($\eta \sim 3$)



Cristián Peña

Caltech






Figure 9. Left: The longitudinal transmittance spectra are shown as a function of wavelength in an expanded scale together with the photo-luminescence spectra for four LSO and LYSO samples before and after the irradiation with integrated doses of 10², 10⁴ and 10⁶ rad. Right: The normalized light output is shown as a function of the integration dose for four long LSO and LYSO samples with PMT (top) and APD (bottom) as the readout devices.






Figure 9. Measurement of the scintillation efficiency of small 1-inch by 1-cm thick samples of scintillator materials exposed to 1Mrad of ^{60}Co irradiation. What is plotted is the ratio of light yield after irradiation to before irradiation for samples as a function of their maximum emission wavelength. The results indicate that for these very small samples, recovery to 90% of original efficiency occurs for essentially all materials. As such these measurements are not particularly informative. More definitive studies are required of actual detector materials (tiles and fibers in needed lengths) in order to draw meaningful conclusions.



Cristián Peña
Caltech

